



**GAMMON TECHNICAL PRODUCTS, INC.**  
 P.O. BOX 400 - 2300 HWY 34  
 MANASQUAN, N.J. 08736

PHONE 732-223-4600  
 FAX 732-223-5778  
 WEBSITE [www.gammontech.com](http://www.gammontech.com)  
 STORE [www.gammontechstore.com](http://www.gammontechstore.com)

**MATCHED-  
 WEIGHT AND  
 COLOR RATING  
 MEMBRANES**

**BULLETIN 73  
 (09-20)**

## COLOR RATING MEMBRANES & MATCHED-WEIGHT MEMBRANES

### APPROVED FOR TESTING JET FUEL PER ASTM METHODS D2276/IP216 AND D5452

Our 37mm, 0.8µm (micrometer) membranes are produced in exact accordance with ASTM Research Report D1012, under which they have been granted approval for use in methods D2276/IP216 as well as Ip216. We manufacture our own plastic monitors using clear Tenite. Our monitors are made to fit all recognized brands of holders including the Gammon MiniMonitor® Kit and the Millipore® Corp. XX64-037-30.


Refill membranes are packaged with a blue paper separator between each one to aid in separation.

Matched-weight refill membranes are packaged with each pair in a glassine envelope. Each pair is matched in weight within 100µg (micrograms) in a closely controlled atmosphere.

Support pads are 34mm in diameter so as to fit into the recess of the plastic monitor beneath the membrane.

Matched-weight monitors have yellow plastic plugs in the inlet port and a red plug in the outlet. Single membrane monitors, used for color rating tests, have blue plugs in the inlet and red plugs in the outlet.

TYPE	MODEL NUMBER	INTERCHANGEABLE WITH MILLIPORE P/N	DESCRIPTION (ALL MEMBRANES 0.8 MICRON)
SINGLE MEMBRANE MONITORS	GTP-1985	MAWP-037-PO	48 PLASTIC MONITORS, EACH WITH A MEMBRANE AND SUPPORT PAD
MATCHED-WEIGHT MONITORS	GTP-1986	MAWP-037-PM	48 PLASTIC MONITORS, EACH WITH A PAIR OF MEMBRANES AND A SUPPORT PAD
EMPTY MONITORS	GTP-1987	NONE	48 EMPTY PLASTIC MONITORS
SINGLE MEMBRANE REFILLS	GTP-1983	AAWP-037-PO	100 MEMBRANES AND 100 SUPPORT PADS
MATCHED-WEIGHT REFILLS	GTP-1984	AAWP-037-PM	50 MATCHED-WEIGHT PAIRS OF MEMBRANES AND 50 SUPPORT PADS



See GamGram No. 25 for filter membrane testing and rating.

MiniMonitor Kit is a registered trademark of Gammon Technical Products. Millipore is a registered trademark of the Millipore Corp.

## MEMBRANE FILING NOTEBOOK AND PAGES

Black 3-ring leatherette notebook holds membrane filling pages. Pages made of heavy-duty plastic, each with space to hold 15 mini-envelopes. Filing pages come in a packet of ten (10) sheets, which is the minimum order.

GTP-1400      Notebook  
GTP-1331      Filing pages



---

## MINI-ENVELOPES

These mini-envelopes are zip-lock clear polyfilm bags made to ensure protection and preservation of filter membranes. They feature a specially prepared white panel for writing data regarding the test such as date, rating, location, etc. Mini-envelopes come in a package of five hundred.

GTP-1267      Mini-envelopes



---

## TWEEZERS

Our tweezers have blunt points and a slip-lock device. They grip and hold membranes without needing to be held tight in your hand.

GTP-2099      Tweezers



---

## 47mm MEMBRANES

We have these filter membranes in stock. They are intended for testing jet fuel for particulate content in accordance with ASTM Method D5452. The membranes are white and will fit all brands and models of filter holders that accept 47mm membranes.

TL-3777      Package of 100 membranes, 0.65µm  
TL-3777-1      Package of 100 membranes, 0.8µm